


Halo Peptide ES-C18

Protein and Peptide Separations



Column: 4.6 x 100 mm; Flow rate: 1.5 mL/min; Temperature: 30° C
A: 0.1% TFA/10% ACN, B: 0.1% TFA/70% ACN
Gradient: 0% to 50% B in 15 min.; Injection volume: 5 μ L



Sample 1

Gly-Tyr, Val-Tyr-Val, Met-enk, Angiotensin II, Leu-enk

Ribonuclease, Porcine Insulin


Halo Peptide ES-C18

Protein and Peptide Separations

Column: 4.6 x 100 mm; Flow rate: 1.5 mL/min; Temperature: 30° C

A: 0.1% TFA/10% ACN, B: 0.1% TFA/70% ACN

Gradient: 15% to 50% B in 15 min.; Injection volume: 5 μ L



Sample 2

Leu-enk

Bovine Insulin, Human Insulin, Cytochrome C, Lysozyme