

Periodic Table of The Analytical Fittings

PORT THREAD	3750 PSI 1/16 M-644-03 / M-650 2.3 In-lbs Headless Micro Knurl PEEK	15000 PSI 1/32 UH-904 3.4 In-lbs Headless Knurl PEEK Blend	2500 PSI 1/8 P-331 / P-350 3-4 In-lbs Standard Knurl PEEK FITTINGS 1 Pressure Rating 2 Tubing Size 3 Fitting / Ferrule Part Number 4 Torque Required 5 Head Style 6 Fitting / Ferrule Material	TOOL N-290 3-4 In-lbs Micro Headless Knurl TOOLS 7 Tool Category 8 Tool Part Number 9 Torque Required 10 Head Style	ICONS Tool Required	For more information, please visit www.biotech.se , or contact your local Biotech office:	30000 PSI 1/32 VHP-700 20 In-lbs 4 mm Hex SS	PORT THREAD										
6-40						Global Headquarters Biotech AB www.biotech.se info@biotech.se +46 (0) 300 569180	North America Biotech USA LLC www.biotech.se bmathieu@biotechusa.us +1 (0) 612 703 5718	Japan BioNik Inc www.bionikinc.com info@bionikinc.com +81 (0) 545 38 9125										
6-32	3750 PSI 1/16 M-660 / M-650 2.3 In-lbs Headless Micro Knurl PEEK	4000 PSI .025 F-125 3.4 In-lbs Standard Micro Knurl PEEK							25000 PSI 1/32 VHP-920 3 In-lbs 4 mm Hex SS / PEEK Blend	30000 PSI 1/32 VHP-900 / VHP-900-01 20 In-lbs 4 mm Hex SS								
10-32	4000 PSI 1/16 F-100 3.4 In-lbs Diamond Knurl PCTFE	5000 PSI 1/16 F-130 3.4 In-lbs Standard Knurl PEEK	5000 PSI 1/32 F-354 / LT-132 4 In-lbs Flush Nut SS / PEEK	5000 PSI 1/32 LT-110 / LT-132 4 In-lbs Standard Knurl PEEK	6000 PSI 190 µm F-200 / F-148 3.4 In-lbs Single Wing Delrin / PCTFE	6000 PSI 360 µm F-200 / F-151 3.4 In-lbs Single Wing Delrin / PCTFE	6000 PSI 1/32 F-200 / F-113 3.4 In-lbs Single Wing Delrin / PEEK	6000 PSI 1/16 F-200 / F-142 3.4 In-lbs Single Wing Delrin / PEEK	9000 PSI 1/16 F-195 / F-192 4 In-lbs Standard Knurl PEEK	9000 PSI 1/16 F-198 / F-192 4 In-lbs Headless Knurl PEEK	9000 PSI 1/16 F-287 / F-192 4 In-lbs Standard Knurl PEEK / PEEK Blend	11000 PSI 1/16 VHP-3200 3.4 In-lbs Diamond Knurl SS / PEEK Blend	12000 PSI 1/16 PK-120BLK 8 In-lbs Standard Knurl PEEK Blend	20000 PSI 1/16 U-400 / U-401 20 In-lbs 1/4" Hex SS	20000 PSI 1/16 U-420 / U-401 20 In-lbs 1/4" Hex SS	25000 PSI 1/16 VHP-320 10 In-lbs 1/4" Hex SS / PEEK Blend	30000 PSI 1/16 VHP-200 / VHP-200-01 20 In-lbs 1/4" Hex SS	
10-32	5000 PSI 1/16 F-120 3.4 In-lbs Standard Knurl PEEK	5000 PSI 1/16 F-127 3.4 In-lbs Standard Knurl PEEK	5000 PSI 1/16 F-354 / LT-100 4 In-lbs Flush Nut SS / PEEK	5000 PSI 1/16 LT-110 / LT-100 4 In-lbs Standard Knurl PEEK	6000 PSI 190 µm F-300 / F-148 3.4 In-lbs Double Wing PEEK / PCTFE	6000 PSI 360 µm F-300 / F-151 3.4 In-lbs Double Wing PEEK / PCTFE	6000 PSI 1/32 F-300 / F-113 3.4 In-lbs Double Wing PEEK	6000 PSI 1/16 F-300 / F-142 3.4 In-lbs Double Wing PEEK	9000 PSI 1/16 F-196 / F-192 4 In-lbs Standard Knurl PEEK	9000 PSI 1/16 F-350 / F-192 4 In-lbs Flush Nut SS / PEEK Blend	9000 PSI 1/16 F-284 / F-192 4 In-lbs Headless Knurl PEEK / PEEK Blend	10000 PSI 1/16 LT-110 / LT-135 7 In-lbs Standard Knurl PEEK / PEEK Blend	16500 PSI 1/32 PK-110 / PK-132 8 In-lbs Standard Knurl PEEK Blend	16500 PSI 1/16 PK-110 / PK-100 8 In-lbs Standard Knurl PEEK Blend	20000 PSI 1/16 U-320 / U-321 20 In-lbs 1/4" Hex SS	25000 PSI 1/16 VHP-325 10 In-lbs 1/4" Hex SS / PEEK Blend	30000 PSI 1/16 VHP-205 / VHP-200-01 20 In-lbs 1/4" Hex SS	
1/4-28	1000 PSI 1/16 P-249 2.3 In-lbs Standard Knurl PEEK	1000 PSI 1/8 P-349 2.3 In-lbs Standard Knurl PEEK	1000 PSI 1/8 P-331 / P-359 3.4 In-lbs Standard Knurl PEEK / ETFE	1350 PSI 1/16 F-356 / P-259 3.4 In-lbs Standard Knurl SS / ETFE	1350 PSI 1/16 P-255 / P-259 3.4 In-lbs Standard Knurl PEEK / ETFE	1500 PSI 1/8 F-156 / P-360 3.4 In-lbs Female Knurl PEEK	1500 PSI 1/8 P-331 / P-360 3.4 In-lbs Standard Knurl PEEK	1850 PSI 1/16 F-356 / P-260 3.4 In-lbs Flush Nut SS / PEEK	1850 PSI 1/16 P-255 / P-260 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/32 F-356 / P-248 3.4 In-lbs Flush Nut SS / ETFE	2500 PSI 1/32 P-255 / P-248 3.4 In-lbs Standard Knurl PEEK / ETFE	2500 PSI 1/16 F-356 / P-250 3.4 In-lbs Flush Nut SS / PEEK	2500 PSI 1/16 P-255 / P-250 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/8 F-156 / P-350 3.4 In-lbs Female Knurl PEEK	2500 PSI 1/8 P-331 / P-350 3.4 In-lbs Standard Knurl PEEK	4500 PSI 1/8 F-364 / LT-200 4 In-lbs Flush Nut SS / PEEK	20000 PSI 1/8 C-235 / C-236 40 In-lbs 5/16" Hex SS / PEEK	
1/4-28	1000 PSI 1/8 P-100 3.4 In-lbs Diamond Knurl PCTFE	1000 PSI 1/8 F-364 / P-359 3.4 In-lbs Flush Nut SS / ETFE	1000 PSI 1/8 P-336 / P-359 3.4 In-lbs Headless Knurl PEEK / ETFE	1350 PSI 1/16 LT-115 / P-259 3.4 In-lbs Standard Knurl PEEK / ETFE	1350 PSI 1/16 P-287 / P-259 3.4 In-lbs Headless Knurl PPS / ETFE	1500 PSI 1/8 P-381 / P-360 3.4 In-lbs Standard Knurl PPS / PEEK	1500 PSI 1/8 F-364 / P-360 3.4 In-lbs Flush Nut SS / PEEK	1500 PSI 1/8 P-336 / P-360 3.4 In-lbs Headless Knurl PEEK	1850 PSI 1/16 LT-115 / P-260 3.4 In-lbs Standard Knurl PEEK	1850 PSI 1/16 P-287 / P-260 3.4 In-lbs Headless Knurl PPS / PEEK	2500 PSI 1/32 LT-115 / P-248 3.4 In-lbs Standard Knurl PEEK / ETFE	2500 PSI 1/32 P-287 / P-248 3.4 In-lbs Headless Knurl PPS / ETFE	2500 PSI 1/16 LT-115 / P-250 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/16 P-287 / P-250 3.4 In-lbs Headless Knurl PPS / PEEK	2500 PSI 1/8 F-364 / P-350 3.4 In-lbs Flush Nut SS / PEEK	2500 PSI 1/8 P-336 / P-350 3.4 In-lbs Standard Knurl PEEK	4500 PSI 1/8 P-336 / P-350 3.4 In-lbs Standard Knurl PEEK	4500 PSI 1/8 C-235 / LT-200 40 In-lbs 5/16" Hex SS / PEEK
1/4-28	1000 PSI 1/8 F-156 / P-359 3.4 In-lbs Female Knurl PEEK / ETFE	1000 PSI 1/8 LT-215 / P-359 3.4 In-lbs Standard Knurl PEEK / ETFE	1000 PSI 1/8 P-387 / P-359 3.4 In-lbs Headless Knurl PPS / ETFE	1350 PSI 1/16 P-232 / P-259 3.4 In-lbs Headless Knurl PEEK / ETFE	1350 PSI 1/16 P-420 / P-259 3.4 In-lbs Standard Knurl PEEK	1500 PSI 1/8 LT-215 / P-360 3.4 In-lbs Standard Knurl PEEK	1500 PSI 1/8 P-387 / P-360 3.4 In-lbs Headless Knurl PPS / PEEK	1850 PSI 1/16 P-232 / P-260 3.4 In-lbs Headless Knurl PEEK	1850 PSI 1/16 P-420 / P-260 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/32 P-232 / P-248 3.4 In-lbs Headless Knurl PEEK / ETFE	2500 PSI 1/32 P-420 / P-248 3.4 In-lbs Female Knurl PEEK / ETFE	2500 PSI 1/16 P-232 / P-250 3.4 In-lbs Headless Knurl PEEK	2500 PSI 1/16 P-420 / P-250 3.4 In-lbs Female Knurl PEEK	2500 PSI 1/8 LT-215 / P-350 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/8 P-387 / P-350 3.4 In-lbs Standard Knurl PPS / PEEK	4500 PSI 1/8 LT-215 / LT-200 4 In-lbs Standard Knurl PEEK	4500 PSI 1/8 LT-210 / LT-200 4 In-lbs Double Wing PEEK	
M6	1000 PSI 1/16 P-229 2.3 In-lbs Standard Knurl PEEK	1000 PSI 1/8 P-329 2.3 In-lbs Standard Knurl PEEK	1000 PSI 2.5 mm P-317 / P-366 3.4 In-lbs Standard Knurl PEEK	1350 PSI 1/16 P-213 / P-259 3.4 In-lbs Headless Knurl PEEK / ETFE	1350 PSI 1/16 P-219 / P-259 3.4 In-lbs Standard Knurl PEEK / ETFE	1350 PSI 1/16 P-217 / P-259 3.4 In-lbs Standard Knurl PEEK / ETFE	2500 PSI 1/32 P-213 / P-248 3.4 In-lbs Headless Knurl PEEK / ETFE	2500 PSI 1/32 P-219 / P-248 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/16 P-213 / P-250 3.4 In-lbs Headless Knurl PEEK	2500 PSI 1/8 P-219 / P-250 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/8 P-317 / P-352 3.4 In-lbs Standard Knurl PPS / PEEK	2500 PSI 1/8 P-319 / P-352 3.4 In-lbs Standard Knurl PEEK	2500 PSI 1/8 P-317 / P-355 3.4 In-lbs Standard Knurl PPS / PCTFE	2500 PSI 1.8 mm P-319 / P-355 3.4 In-lbs Standard Knurl PEEK / PCTFE	2500 PSI 1.8 mm P-319 / P-355 3.4 In-lbs Standard Knurl PEEK / PCTFE	5000 PSI 2.0 mm P-357 / P-357-2 3.4 In-lbs Standard Knurl PEEK	9000 PSI 1/16 F-197 / F-192 4 In-lbs Standard Knurl PEEK	20000 PSI 1/16 U-450 / U-401 20 In-lbs 5/16" Hex SS

Copyright © 2015 IDEX Health & Science LLC. IDEX Health & Science LLC is a Unit of IDIX Corporation. | IDI015

HAND TIGHTENING/EXTENDER TOOLS

F-345 3.4 In-lbs Flush Nut-10-32	F-346 3.4 In-lbs Flush Nut 1/4-28	N-290 3.4 In-lbs Micro Headless Knurl	P-277 3.4 In-lbs Standard Micro Knurl	P-297 3.4 In-lbs Headless Knurl	P-299 3.4 In-lbs Standard Knurl	P-399 3.4 In-lbs Standard Micro Knurl	VHP-1001 10 In-lbs 1/4" Hex	VHP-9000 2-10 In-lbs 4 mm Hex	VHP-1000 10 In-lbs 1/4" Hex	VHP-2000 14 In-lbs 1/4" Hex	VHP-3000 10 In-lbs 8 mm Hex	P-1000 4 In-lbs Standard Knurl
-------------------------------------	--------------------------------------	--	--	------------------------------------	------------------------------------	--	--------------------------------	----------------------------------	--------------------------------	--------------------------------	--------------------------------	-----------------------------------

ADJUSTABLE TORQUE DRIVER AND TORQUE DRIVER EXTENDER TOOLS

VHP-4000 3-10 In-lbs Double Drive	N-291 2-10 In-lbs Micro Headless Knurl	P-268 2-10 In-lbs 1/4" Hex	P-278 2-10 In-lbs Female Nut Knurl	P-279 2-10 In-lbs Micro Knurl	P-291 2-10 In-lbs Standard Knurl	P-292 2-10 In-lbs Headless Knurl	F-347 2-10 In-lbs Flush Nut 1/4-28
--------------------------------------	---	-------------------------------	---------------------------------------	----------------------------------	-------------------------------------	-------------------------------------	---------------------------------------