

DEGASi® PLUS MICRO

Are you working with very low flow rates? Do you want to minimize the dead volumes inside your system? Then DEGASi® MICRO is your number one choice of degasser. With only 100 µl internal volume you minimize the internal volume while still keeping excellent degassing efficiency up to approximately 0.5 ml/min.

"an outstanding degasser

for very low flow rates..."

PRODUCT FEATURES

- Ultra-high degassing efficiency
- Fits every system
- Long life expectancy 5+ years continuous (24/7) running capacity • Continuous vacuum system
- Easy to prime
- Almost silent
- Closed loop vacuum control means constant vacuum (variable • Advanced error and leak checking RPM)
- ZHCR® patented control eliminates CE & RoHS compliant baseline fluctuations

	.	
Part Number	Number of Channels	Internal Volume
0003-6351-S	1	100 μΙ
0003-6352-S	2	100 μΙ
0003-6353-S	3	100 μΙ
0003-6354-S	4	100 μΙ
0003-6355-S	5	100 μΙ
0003-6356-5	6	100 μΙ

- Excellent chemical compatibility flow path
- Biocompatible flow path
- monitoring to ensure that optimal operational conditions are maintained
- functions



GENERAL SPECIFICATIONS

Degassing Channel Tubing:

Systec AF™ (0.045" ID)

Maximum Channel Pressure:

70 PSI

Wetted Materials:

Systec AF™, PEEK, Glass-filled PTFE

Liquid Connection:

1/4"-28 UNF threaded flat-bottom port

Size (L x H x W):

245 x 111 x 59 mm

POWER REQUIREMENTS **USING SUPPLIED AC ADAPTER**

24 V Power supply included. Required input to power supply: 110-240 V 50/60 HZ

Wall Sockets: 4 supplied with AC Adapter, Interchangeable: North America, Japan, U.K., Continental Europe, Australia

16